



Anti-Eotaxin (aa 83-97) polyclonal antibody (DPAB-DC4207)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 83-97 of human Eotaxin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate),
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The eotaxins are a CC chemokine subfamily of eosinophil chemotactic proteins.
Keywords	Eotaxin

GENE INFORMATION

Entrez Gene ID [6356](#)

UniProt ID [P51671](#)
