



Anti-Alligator ER-alpha polyclonal antibody (DPAB-DC4213)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to alligator Estrogen receptor-alpha.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Crocodile
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	25 µg
Buffer	In PBS (0.1% proclin, 2.0% Block Ace)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Estrogen receptor alpha (ER-α), also known as NR3A1 (nuclear receptor subfamily 3, group A, member 1), is one of two main types of estrogen receptor, a nuclear receptor that is activated by the sex hormone estrogen. In humans, ER-α is encoded by the gene
Keywords	ER-α; ER-alpha