



Anti-E. coli O26:B6, O26:B6 Polyclonal antibody (DPAB-DC4219)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	E. coli O26:B6, O26:B6
Immunogen	Escherichia coli O26:B6 total proteins.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction Escherichia coli (E. coli) is Gram negative enteric bacteria that live in the human gut. Its

presence in nature is an indication of human faecal pollution. E. coli also causes human disease such as urinary tract infections and neonatal meningitis. Infecti

Keywords

E. coli O26:B6; Escherichia coli O26:B6
