



# Anti-E. coli Polyclonal antibody (DPAB-DC4221)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Reactive to all E. coli species in Immunogen, plus O157-H7 tested by immunoblot. Reacts with more than 50 proteins by Western blot.
<b>Target</b>	E. coli
<b>Immunogen</b>	Native Escherichia coli: TOP 10F, HB101, JM109, BL21, DH5 alpha and K12.
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	E. coli
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA, Dot, EIA,
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	In citrate buffer
<b>Preservative</b>	None
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Escherichia coli (commonly abbreviated E. coli) is a Gram-negative, facultatively anaerobic, rod-shaped bacterium of the genus Escherichia that is commonly found in the lower intestine of
---------------------	---

warm-blooded organisms (endotherms). Most E. coli strains are harml

---

**Keywords**

E. coli; Escherichia coli

---