



Anti-E. coli Polyclonal antibody (DPAB-DC4222)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reactive to all E. coli species in Immunogen, plus O157-H7 tested by immunoblot. Reacts with more than 50 proteins by Western blot.
Target	E. coli
Immunogen	Native Escherichia coli: TOP 10F, HB101, JM109, BL21, DH5 alpha and K12.
Source/Host	Goat
Species Reactivity	E. coli
Conjugate	Unconjugated
Applications	WB, ELISA, Dot, EIA,
Format	Liquid
Size	1 ml
Buffer	In citrate-phosphate buffer, pH 7.0-7.4
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Escherichia coli (commonly abbreviated E. coli) is a Gram-negative, facultatively anaerobic, rod-shaped bacterium of the genus Escherichia that is commonly found in the lower intestine of
---------------------	---

warm-blooded organisms (endotherms). Most E. coli strains are harml

Keywords

E. coli; Escherichia coli
