



# Anti-Ethinylestradiol polyclonal antibody (DPAB-DC4225)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	1, 3, 5(10)-ESTRATRIEN-17a-ETHYNYL-3, 17b-DIOL-6-ONE 6-CARBOXYMETHYLOXIME.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Ethinyl estradiol (EE), also sometimes written as ethinylestradiol, ethynyl estradiol, or ethinyl estradiol, is a derivative of 17β-estradiol (E2), the major endogenous estrogen in humans. EE is an orally bioactive estrogen used in many formulations of co
Keywords	EE; Ethinylestradiol