



Anti-Enterovirus VP1 Polyclonal antibody (DPAB-DC4227)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Enterovirus VP1
Immunogen	EV71 VP1 (1 a.a. ~ 297 a.a) full-length human protein.
Source/Host	Mouse
Species Reactivity	Enterovirus
Conjugate	Unconjugated
Applications	WB (Transfected lysate), Detection Antibody,
Format	Liquid
Size	50 µg
Buffer	In PBS, pH 7.4
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	EV71 VP1; EV71 VP1 purified MaxPab mouse
----------	--