



Anti-Factor XIII polyclonal antibody (DPAB-DC4229)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody recognizes both alpha and beta subunits of Factor XIII with minimal cross reactivity with other human serum proteins.
Immunogen	Purified human Factor XIII.
Source/Host	Sheep
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ICC, IF, ELISA, ID, Immunoelectron Microscopy,
Format	Liquid
Size	1 ml
Buffer	In serum
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords Factor XIII

GENE INFORMATION

Entrez Gene ID [14060](#)

UniProt ID [Q07968](#)
