



Anti-FANCC (aa 550-622) polyclonal antibody (DPAB-DC4231)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to amino acids 550-622 of FANCC.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Transfected lysate),
Format	Liquid
Size	100 µl
Buffer	In antiserum
Preservative	None
Storage	Store at -20°C or -80°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	FANCC; FANCG
----------	--------------

GENE INFORMATION

Entrez Gene ID [2176](#)

UniProt ID [Q00597](#)
