



# Anti-FBP1 (native protein) polyclonal antibody (DPAB-DC4233)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Native purified bovine folate binding protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine
Conjugate	Unconjugated
Applications	WB, IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	FBP1
----------	------

## GENE INFORMATION

Entrez Gene ID [2203](#)

---

UniProt ID [P09467](#)

---