



Anti-FBX4 (N-terminal) polyclonal antibody (DPAB-DC4234)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antiserum is directed against human and mouse FBX4.
Immunogen	A synthetic peptide corresponding to N-terminus of human and mouse FBX4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB (Cell lysate), IF, ELISA,
Format	Liquid
Size	100 µl
Buffer	In 20 mM KH ₂ PO ₄ , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	FBX4
----------	------

GENE INFORMATION

Entrez Gene ID [26272](#)

UniProt ID [B3KNA0](#)
