



Anti-FeLV P27 Polyclonal antibody (DPAB-DC4235)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	FeLV p27 (Western blot). Non-reactive with normal feline serum.
Target	FeLV P27
Immunogen	Native purified Feline Leukemia Virus p27.
Source/Host	Goat
Species Reactivity	FeLV
Conjugate	Unconjugated
Applications	WB, IF,
Format	Liquid
Size	1 ml
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Feline leukemia virus (FeLV) is a retrovirus that infects cats. FeLV can be transmitted from infected cats when the transfer of saliva or nasal secretions is involved. If not defeated by the animal's immune system, the virus can cause diseases which can b
--------------	---

Keywords

FeLV p27; Feline Leukemia Virus p27
