



Anti-Ferritin (native) polyclonal antibody (DPAB-DC4238)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified ferritin from rat liver.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (0.02% sodium azide, 50% glycerol (v/v))
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Ferritin is a ubiquitous intracellular protein that stores iron and releases it in a controlled fashion. The protein is produced by almost all living organisms, including algae, bacteria, higher plants, and animals. In humans, it acts as a buffer against
---------------------	--

Keywords

Ferritin

GENE INFORMATION

Entrez Gene ID

[2495](#)

UniProt ID

[P02794](#)
