



Anti-Fibrinogen fragment D (native) polyclonal antibody (DPAB-DC4240)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with intact fibrinogen and fibrinogenin from plasma. No cross-reaction with human fibrinogen fragment E (when assessed with double diffusion/electrophoresis). No cross reactivity with human serum.
Immunogen	Native purified plasmin digestion product D (high molecular weight) of human fibrinogen.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, ICC, IP, ELISA,
Format	Liquid
Size	1 ml
Buffer	In antiserum
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords Fibrinogen fragment D