



# Anti-Fibrinogen fragment E (native) polyclonal antibody (DPAB-DC4241)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Specificity** This antibody reacts with intact fibrinogen and fibrinogenin from human plasma. No cross-reactivity with human fibrinogen fragment D (when assessed with double diffusion/electrophoresis).

**Immunogen** Native purified plasmin digestion product E of human fibrinogen.

**Isotype** IgG

**Source/Host** Rabbit

**Species Reactivity** Human

**Conjugate** Unconjugated

**Applications** WB, ICC, IP, ELISA,

**Format** Liquid

**Size** 1 ml

**Buffer** In antiserum

**Preservative** None

**Storage** Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

**Keywords** Fibrinogen fragment E