



Anti-Fibrinogen (native) polyclonal antibody (DPAB-DC4246)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified mouse fibrinogen.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB, IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Fibrinogen (factor I) is a glycoprotein in vertebrates that helps in the formation of blood clots. It consists of a linear array of 3 nodules held together by a very thin thread which is estimated to have a diameter of from 8 to 15 Å. The two end nodules
Keywords	Fibrinogen

GENE INFORMATION

Entrez Gene ID [2243](#)

UniProt ID [P02671](#)
