



Anti-Fibrinopeptide B (C-terminal) polyclonal antibody (DPAB-DC4247)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Glp-Gly-Val-Asn-Asp-Asn-Glu-Glu-Gly-Phe-Phe-Ser-Ala-Arg conjugate to KLH via carboxyl group.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purification
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (0.02% sodium azide, 50% glycerol (v/v))
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	Fibrinopeptide B
----------	------------------

GENE INFORMATION

Entrez Gene ID [2244](#)

UniProt ID [P02675](#)
