



Anti-Fibronectin (native) polyclonal antibody (DPAB-DC4248)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody has a single antigen reactivity with minimal cross reactivity with other human plasma proteins.
Immunogen	Native purified from human plasma Fibronectin.
Source/Host	Sheep
Species Reactivity	Human, Rat
Conjugate	Unconjugated
Applications	ICC, IF, ELISA, RIA, ID, Immunoelectron Microscopy,
Format	Liquid
Size	25 µg, 100 µg
Buffer	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. Reconstitute at 0.2 mg/mL in sterile PBS.
Preservative	None
Storage	Store at 4°C.

BACKGROUND

Introduction	Fibronectin is a high-molecular weight (~440kDa) glycoprotein of the extracellular matrix that binds to membrane-spanning receptor proteins called integrins. Similar to integrins, fibronectin binds extracellular matrix components such as collagen, fibrin,
---------------------	---

Keywords

Fibronectin

GENE INFORMATION

Entrez Gene ID

[2335](#)

UniProt ID

[P02751](#)
