



# Anti-FT (native) polyclonal antibody (DPAB-DC4258)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	By immunoelectrophoresis, This antibody reacts with human ferritin. There is no reaction with normal human serum.
<b>Immunogen</b>	Native purified human FT.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	1 ml
<b>Buffer</b>	In PBS, pH 7.2 (0.1% sodium azide)
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Keywords</b>	FT
-----------------	----