



# Anti-GABBR2 (N-terminal) polyclonal antibody (DPAB-DC4260)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Specificity	Reacts to the N-terminal region of the GABA BR2 protein.
Immunogen	A synthetic peptide corresponding to N-terminus of GABBR2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µl
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for short term. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	GABBR2
----------	--------

# GENE INFORMATION

Entrez Gene ID [9568](#)

---

UniProt ID [H9NIL8](#)

---