



# Anti-E. coli GADE Polyclonal antibody (DPAB-DC4261)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	E. coli GADE
Immunogen	A synthetic peptide corresponding to N-terminus of Escherichia coli gadE.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	GADE
----------	------