



# Anti-Galantide polyclonal antibody (DPAB-DC4262)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide (conjugated with KLH) corresponding to galantide. The sequence is GWTLNSAGYLLGPQQFFGLM
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IP, EIA,
<b>Format</b>	Liquid
<b>Size</b>	150 µg
<b>Buffer</b>	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Galantide is a galanin antagonist consisting of aa 1-13 of galanin and the C-terminal fragment of bradykinin. Galantide reversibly blocks the neuronal activities of galanin. Galanin has different activities such as contraction of smooth muscle of the gast
<b>Keywords</b>	Galantide