



# Anti-Gastrin 1 (C-terminal) polyclonal antibody (DPAB-DC4263)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Glp-Gly-Pro-Trp-Leu-Glu-Glu-Glu-Glu-Glu-Ala-Tyr-Gly-Trp-Met-Asp-Phe-NH <sub>2</sub> conjugate to KLH via carboxyl group.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purification
Conjugate	Unconjugated
Applications	IP, EIA,
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (0.02% sodium azide, 50% glycerol (v/v))
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	This peptide, mainly expressed in the gut, activates gastric acid secretion and stimulates lectin-like protein Reg expression via activation of PKC and RhoA. It is used as a marker for gastritis.
--------------	---

**Keywords**Gastrin-1; Gastrin 1

---

**GENE INFORMATION****Entrez Gene ID**[2520](#)

---

**UniProt ID**[P01350](#)

---