



Anti-Gastrin Releasing Peptide polyclonal antibody (DPAB-DC4264)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with Tg) corresponding to frog Gastrin releasing peptide.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Frog, Human, Mouse, Rat
Conjugate	Unconjugated
Applications	IHC-P,
Format	Liquid
Size	250 µl
Buffer	In PBS, pH 7.3 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	In humans, gastrin is a peptide hormone that stimulates secretion of gastric acid (HCl) by the parietal cells of the stomach and aids in gastric motility. It is released by G cells in the antrum of the stomach (the portion of the stomach adjacent the pylo
Keywords	Gastrin Releasing Peptide

GENE INFORMATION

Entrez Gene ID [2922](#)

UniProt ID [P07492](#)
