



Anti-G. lamblia Polyclonal antibody (DPAB-DC4275)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Giardia lamblia intact cysts. Cross-reacts with a rod-shaped microorganism found in rodent feces.
Target	G. lamblia
Immunogen	Native purified from Giardia lamblia intact cysts.
Source/Host	Goat
Species Reactivity	G. lamblia
Conjugate	Unconjugated
Applications	IF, ELISA, EIA,
Format	Liquid
Size	1 ml
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Giardia lamblia is a flagellated protozoan parasite that colonizes and reproduces in the small intestine, causing giardiasis. The parasite attaches to the epithelium by a ventral adhesive disc,
---------------------	---

and reproduces via binary fission. Giardiasis does not spread

Keywords

G. lamblia; Giardia lamblia
