



# Anti- Gibberellic acid Polyclonal antibody (DPAB-DC4276)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	Gibberellic acid
Immunogen	Gibberellic acid conjugated with carrier proteins.
Source/Host	Rat
Species Reactivity	N/A
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	IHC, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Gibberellic acid (also called Gibberellin A3, GA, and GA3) is a hormone found in plants. Its chemical formula is C <sub>19</sub> H <sub>22</sub> O <sub>6</sub> . When purified, it is a white to pale-yellow solid. Gibberellic acid is a simple gibberellin, a pentacyclic diterpene acid promoting
--------------	---

**Keywords**

GA; Gibberellic acid

---