



Anti-Glucagon (native) polyclonal antibody (DPAB-DC4281)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified Glucagon from pig pancreatic (conjugated with BSA).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Pig, Rat
Conjugate	Unconjugated
Applications	IHC,
Format	Liquid
Size	250 µl
Buffer	In PBS, pH 7.3 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Glucagon is a peptide hormone, produced by alpha cells of the pancreas, that raises the concentration of glucose in the bloodstream. Its effect is opposite that of insulin, which lowers the glucose concentration. The pancreas releases glucagon when the co
Keywords	Glucagon

GENE INFORMATION

Entrez Gene ID [2641](#)

UniProt ID [P01275](#)
