



Anti-R. glucoamylase Glucoamylase 1 Polyclonal antibody (DPAB-DC4282)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Target | R. glucoamylase Glucoamylase 1 |
| Immunogen | Native purified Rhizopus Glucoamylase 1. |
| Isotype | IgG |
| Source/Host | Goat |
| Species Reactivity | R. glucoamylase |
| Conjugate | Unconjugated |
| Applications | WB (Recombinant protein), ELISA, |
| Format | Lyophilized |
| Size | 100 µg |
| Buffer | Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide) |
| Preservative | 0.01% Sodium Azide |
| Storage | Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing. |

BACKGROUND

| | |
|---------------------|---|
| Introduction | Glucan 1,4-alpha-glucosidase (EC 3.2.1.3, glucoamylase, amyloglucosidase, gamma-amylase, lysosomal alpha-glucosidase, acid maltase, exo-1,4-alpha-glucosidase, glucose amylase, |
|---------------------|---|

gamma-1,4-glucan glucohydrolase, acid maltase, 1,4-alpha-D-glucan glucohydrolas

Keywords

Glucoamylase 1
