



# Anti-GM polyclonal antibody (DPAB-DC4290)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to human GM-CSF.
Source/Host	Sheep
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, EIA,
Format	Liquid
Size	1 ml
Buffer	In 20 mM phosphate buffer, 150 mM NaCl, pH 7.2 (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	GM
----------	----