



Anti-GNB5 polyclonal antibody (DPAB-DC4291)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to the D1 blade of GNB5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB (Recombinant protein),
Format	Liquid
Size	100 µg
Buffer	In PBS (0.08% sodium azide)
Preservative	0.08% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	GNB5
----------	------

GENE INFORMATION

Entrez Gene ID	10681
UniProt ID	O14775