



Anti-GnRHa polyclonal antibody (DPAB-DC4292)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	GnRHa (des-Gly10, D-Ala6, Pro9-NHEt)-LHRH, pGlu-His-Trp-Ser-Tyr-D-Ala-Leu-Arg Pro-NHEt.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA,
Format	Lyophilized
Size	200 µl
Buffer	Lyophilized from serum
Preservative	None
Storage	Store at -20°C on dry atmosphere. After reconstitution with 200 uL deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	GnRHa
-----------------	-------

GENE INFORMATION

Entrez Gene ID [2796](#)

UniProt ID [P01148](#)
