



Anti-Gonadotropin (C-terminal) polyclonal antibody (DPAB-DC4293)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to C-terminus of human Gonadotropin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.1% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Gonadotropins (or glycoprotein hormones) are protein hormones secreted by gonadotrope cells of the anterior pituitary of vertebrates. This is a family of proteins, which include the mammalian hormones follicle-stimulating hormone (FSH), luteinizing hormone
Keywords	Gonadotropin