



Anti-H. alvei Polyclonal antibody (DPAB-DC4303)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	H. alvei
Immunogen	Hafnia alvei total protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	H. alvei
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB, ELISA,
Format	Lyophilized
Size	100 µl
Buffer	Lyophilized from PBS, pH 7.4
Preservative	None
Storage	Store at 4°C on dry atmosphere. After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction Hafnia is the genus of the Enterobacteriaceae family whose only species is the Gram-negative,

facultatively anaerobic, rod-shaped bacterium *Hafnia alvei*.

H. alvei is a commensal of the human gastrointestinal tract and not normally pathogenic, but may cau

Keywords

H. alvei; *Hafnia alvei*
