



# Anti-Hallen (aa 78-91) polyclonal antibody (DPAB-DC4304)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Hallen is a novel BRCT domain containing protein.
<b>Immunogen</b>	A synthetic peptide corresponding to amino acids 78-91 of human Hallen. The sequence is KSGRWLDETTYEWGC
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In 20 mM KH <sub>2</sub> PO <sub>4</sub> , 150 mM NaCl, pH 7.2 (0.01% sodium azide)
<b>Preservative</b>	0.01% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	Hallen is a novel BRCT domain containing protein.
---------------------	---

