



Anti-HAP40 polyclonal antibody (DPAB-DC4305)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody reacts with residues 314-325 [DGHGQDTSGQLP] of the mouse HAP40 protein (co-immunoprecipitates Huntingtin protein).
Immunogen	A synthetic peptide corresponding to mouse HAP40.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for short term. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords HAP40 ; HAP40

GENE INFORMATION

Entrez Gene ID [8263](#)

UniProt ID [P23610](#)
