



# Anti-HB polyclonal antibody (DPAB-DC4306)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Cross-reactivity to other species may exist. Antibody activity to serum proteins is not present.
<b>Immunogen</b>	Human HB.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA, RIA, ID,
<b>Format</b>	Liquid
<b>Size</b>	2 ml
<b>Buffer</b>	In 20 mM sodium phosphate buffer, 140 mM NaCl, pH 7.3 (0.05% sodium azide)
<b>Preservative</b>	0.05% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Keywords</b>	HB
-----------------	----