



Anti-H. pylori Polyclonal antibody (DPAB-DC4310)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	H. pylori
Immunogen	Native purified Helicobacter pylori antigens.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	H. pylori
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	1 mg
Buffer	In 15 mM KPO4 buffer, 0.85% NaCl, pH 7.4 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Helicobacter pylori (H. pylori), previously named Campylobacter pylori, is a Gram-negative, microaerophilic bacterium found in the stomach. It was identified in 1982 by Australian scientists Barry Marshall and Robin Warren, who found that it was present i
--------------	---

Keywords

H. pylori antigen; Helicobacter pylori antigens
