



Anti-H. pylori Polyclonal antibody (DPAB-DC4311)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	H. pylori
Immunogen	Native from Helicobacter pylori.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	H. pylori
Conjugate	Unconjugated
Applications	IHC,
Format	Liquid
Size	100 µg
Buffer	In buffer containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Helicobacter pylori (H. pylori), previously named Campylobacter pylori, is a Gram-negative, microaerophilic bacterium found in the stomach. It was identified in 1982 by Australian scientists Barry Marshall and Robin Warren, who found that it was present i
--------------	---

Keywords

H. pylori; Helicobacter pylori
