



Anti-H. pylori Polyclonal antibody (DPAB-DC4312)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	H. pylori
Immunogen	Purified Helicobacter pylori.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	H. pylori
Purification	Immunoaffinity purification
Conjugate	Unconjugated
Applications	IHC-P,
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	In 10 mM PBS, pH 7.4 (1% BSA, 0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction

Helicobacter pylori (*H. pylori*), previously named *Campylobacter pylori*, is a Gram-negative, microaerophilic bacterium found in the stomach. It was identified in 1982 by Australian scientists Barry Marshall and Robin Warren, who found that it was present i

Keywords

H. pylori; *Helicobacter pylori*
