



Anti-H. pylori Polyclonal antibody (DPAB-DC4314)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	H. pylori
Immunogen	Native purified HPYR1 from Helicobacter pylori strain CH-20 429.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	H. pylori
Conjugate	Unconjugated
Applications	IHC-P,
Format	Lyophilized
Size	1 ml
Buffer	Lyophilized from tissue culture supernatant (15 mM sodium azide)
Preservative	15mM Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with steril distilled water, Store at 4°C for at least two months. For long term storage store at -20°C or below. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Helicobacter pylori (H. pylori), previously named Campylobacter pylori, is a Gram-negative, microaerophilic bacterium found in the stomach. It was identified in 1982 by Australian
---------------------	---

scientists Barry Marshall and Robin Warren, who found that it was present i

Keywords

H. pylori; Helicobacter pylori
