



Anti-IAV H5N1 Hemagglutinin Polyclonal antibody (DPAB-DC4319)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Hemagglutinin protein from the H5N1 strain of avian influenza (Middle region).
Target	IAV H5N1 Hemagglutinin
Immunogen	A synthetic peptide corresponding to internal region 10 amino acids of Influenza A H5N1 (Avian) Hemagglutinin.
Source/Host	Goat
Species Reactivity	IAV
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Influenza A virus subtype H5N1, also known as A(H5N1) or simply H5N1, is a subtype of the influenza A virus which can cause illness in humans and many other animal species. A bird-adapted strain of H5N1, called HPAI A(H5N1) for highly pathogenic avian inf
---------------------	---

Keywords

Hemagglutinin (H5N1)
