



Anti-HBV Core Protein Polyclonal antibody (DPAB-DC4323)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	HBV Core Protein
Immunogen	Native purified virus hepatitis B virus core antigen.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HBV
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr,
Format	Liquid
Size	500 µl
Buffer	In 10 mM phosphate buffer, 150 mM NaCl, pH 7.4 (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Hepatitis B virus, abbreviated HBV, is a species of the genus Orthohepadnavirus, which is likewise a part of the Hepadnaviridae family of viruses. This virus causes the disease hepatitis B.
---------------------	---

Keywords

HBV core antigen; Hepatitis B virus core antigen
