



# Anti-HBV Polyclonal antibody (DPAB-DC4324)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	HBV
Immunogen	A synthetic peptide corresponding to amino acids 84-96 of HBV surface antigen.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HBV
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	Hepatitis B virus, abbreviated HBV, is a species of the genus Orthohepadnavirus, which is likewise a part of the Hepadnaviridae family of viruses. This virus causes the disease hepatitis B.
Keywords	HBV; Hepatitis B virus