



Anti-HBV type ad, ay Surface Antigen Polyclonal antibody (DPAB-DC4325)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	HBsAg specific to ad/ay. Does not react with human serum proteins.
Target	HBV type ad, ay Surface Antigen
Immunogen	Native purified hepatitis B surface antigen (HBsAg) ad/ay.
Source/Host	Horse
Species Reactivity	HBV
Conjugate	Unconjugated
Applications	WB, IHC-Fr, IF, IP, Dot, EIA,
Format	Liquid
Size	1 mg
Buffer	In PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C.

BACKGROUND

Introduction	Hepatitis B virus, abbreviated HBV, is a species of the genus Orthohepadnavirus, which is likewise a part of the Hepadnaviridae family of viruses. This virus causes the disease hepatitis B.
---------------------	---

Keywords

HBV surface antigen; Hepatitis B virus surface antigen (ad/ay)
