



# Anti-HBV type ad, ay Surface Antigen Polyclonal antibody (DPAB-DC4326)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Reacts only with Hepatitis B surface antigen. Non-reactive with normal human serum.
<b>Target</b>	HBV type ad, ay Surface Antigen
<b>Immunogen</b>	Mixture of subtypes ad and ay hepatitis B surface antigen purified from human serum.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	HBV
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

**Introduction** Hepatitis B virus, abbreviated HBV, is a species of the genus Orthohepadnavirus, which is

likewise a part of the Hepadnaviridae family of viruses. This virus causes the disease hepatitis B.

---

**Keywords**

HBV surface antigen; Hepatitis B virus surface antigen (ad/ay)

---