



Anti-HBV Surface Antigen Polyclonal antibody (DPAB-DC4327)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody recognizes human Hepatitis B surface antigen.
Target	HBV Surface Antigen
Immunogen	Native virulates hepatitis B virus surface antigen.
Source/Host	Goat
Species Reactivity	HBV
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Size	100 µg
Buffer	In 10 mM PBS, pH 7.2
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Hepatitis B virus, abbreviated HBV, is a species of the genus Orthohepadnavirus, which is likewise a part of the Hepadnaviridae family of viruses. This virus causes the disease hepatitis B.
---------------------	---

Keywords

HBV surface antigen; Hepatitis B virus surface antigen
