



# Anti-HCV NS5b Polyclonal antibody (DPAB-DC4333)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	HCV NS5b antigen. Recognizes NS5b part of the NS5ab recombinant antigen.
<b>Target</b>	HCV NS5b
<b>Immunogen</b>	Recombinant protein corresponding to hepatitis C virus NS5b.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	HCV
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IP, ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS, pH 7.2 (0.01% sodium azide)
<b>Preservative</b>	0.01% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

**Introduction** Hepatitis C virus (HCV or sometimes HVC) is a small (55–65 nm in size), enveloped, positive-

sense single-stranded RNA virus of the family Flaviviridae. Hepatitis C virus is the cause of hepatitis C in humans.

---

<b>Keywords</b>	HCV NS5b; Hepatitis C virus NS5b
-----------------	----------------------------------

---