



Anti-HCV Polyclonal antibody (DPAB-DC4338)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	NS-5a and b.
Target	HCV
Immunogen	Recombinant protein corresponding to NS-5 a and b (100 % of a and 70 % of b) of hepatitis C virus (genotype 1a).
Source/Host	Goat
Species Reactivity	HCV
Conjugate	Unconjugated
Applications	WB, ELISA, EIA,
Format	Liquid
Size	1 ml
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Hepatitis C virus (HCV or sometimes HVC) is a small (55–69 nm in size), enveloped, positive-sense single-stranded RNA virus of the family Flaviviridae. Hepatitis C virus is the cause of hepatitis C in humans.
---------------------	--

Keywords

HCV; Hepatitis C virus
