



Anti-HGV Polyclonal antibody (DPAB-DC4342)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	HGV
Immunogen	A synthetic peptide corresponding to amino acids 2268-2276 of HGV. The sequence is CDKCEARQE
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HGV
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	GB virus C (GBV-C), formerly known as hepatitis G virus (HGV) and also known as HPgV is a virus in the Flaviviridae family and a member of the Pegivirus genus, is known to infect humans, but is not known to cause human disease. There have been reports tha
---------------------	---

Keywords

HGV; Hepatitis G virus
