



Anti-HSV type 42006 Polyclonal antibody (DPAB-DC4344)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reacts with Herpes simplex virus types 1 and 2, ICPs and late structural (virion) antigens. Does not react with HEp-2 cells.
Target	HSV type 42006
Immunogen	Native from human Herpes Simplex Virus types 1 and 2 Strain F.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HSV
Conjugate	Unconjugated
Applications	IF,
Format	Liquid
Size	100 µl
Buffer	In 10 mM PBS, pH 7.2 (0.1% sodium azide)
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction

Herpes simplex virus 1 and 2 (HSV-1 and HSV-2), also known as human herpesvirus 1 and 2 (HHV-1 and HHV-2), are two members of the herpesvirus family, Herpesviridae, that infect humans. Both HSV-1 (which produces most cold sores) and HSV-2 (which produces

Keywords

HSV-1/2; Herpes simplex virus type 1/2
